

HistoPot

Laboratory chemicals and products-Clinical chemistry and hematology products-Fine & industrial chemicals

Certificate of Analysis

Product Name 10% v/v Buffered Formalin
Formalin Solution 10% (v/v) [=4% (w/v)] neutralized
HISTO GRADE
Fixative for histology

HistoPot Unit 2.5 LTP

Lot No 29081912944

Expiry Date 2022-07

This solution is phosphate buffered at 50 mmol/l.

Test	Specification	Actual value
Formaldehyde (HCHO)	38 - 42 g/l	38 - 42 g/l
pH at 25°C	6.8 - 7.2	6.8 - 7.2
Sodium (Na)	61 - 71 mmol/l	61 - 71 mmol/l

Signature of approval:



Mr. Marcus Nash
Manufacturing Scientist
Serosep Ltd.