

Certificate of Performance

Phosphate 2 (Phos 2) Reagent

Product Code **KR107**
 Lot Number **4318**



Date 17.02.2020


Expiry Date 17.02.2021

| Reagent | Lot Number |
|-----------------|-------------|
| Phos 1 | A19-214 |
| Phos 2 | 4318 |
| Test Method | WC SOP 21 |
| Test Instrument | Konelab 20 |
| Serial Number | P4919028 |

| Standard Value | Lot # | Abs. |
|----------------|-----------|-------|
| 0 | | 0.001 |
| 0.5 | P20 09 19 | 0.149 |
| 1 | P20 09 19 | 0.298 |
| 1.5 | P20 09 19 | 0.458 |
| 2 | P20 09 19 | 0.586 |

*serial dilution

Calibration Curve



| Test | Result | Specification | Pass / Fail |
|-------------|------------------|-----------------|-------------|
| Intercept | 0.003 | -0.005 to 0.010 | Pass |
| Sensitivity | Abs/ppm 0.296 | 0.260 to 0.320 | Pass |
| Regression | 0.999 | 0.990 to 1.001 | Pass |

Notes
 Store between 2 & 8 °C
 *Stable on Board 5 days
 Barcode ID = 107
 *Based on results obtained in our Research Facility

Alan Coonerty
 A. Coonerty

An ISO 9001:2008 Registered Company
 Annacotty Business Park, Annacotty, Limerick, Ireland.
 Telephone 061 358190 Facsimile 061 358191

Hazard Identification

Labelling according to Regulation (EC) No 1272/2008
 This reagent is classed as non hazardous.