

Certificate of Performance

TON 1

Product Code KR103

Date 13.03.2020


Expiry Date 12.03.2021

| Reagent | Lot Number |
|-----------------|-------------|
| TON 1 | 4350 |
| TON 2 | 4362 |
| TON 3 | 4353 |
| Test Method | WC SOP 23 |
| Test Instrument | Konelab 20 |
| Serial Number | P4919028 |

| Std. Value | Lot # | Abs. |
|------------|---------------|-------|
| 0.00 | | 0.001 |
| 2.50 | TON 100 01 20 | 0.170 |
| 5.00 | TON 100 01 20 | 0.327 |
| 7.69 | TON 100 01 20 | 0.495 |
| 10.00 | TON 100 01 20 | 0.619 |

*serial dilution

Calibration Curve



| Test | Result | Specification | Pass / Fail |
|------------------|---------|------------------|-------------|
| Intercept | 0.001 | -0.005 to 0.01 | Pass |
| Sensitivity | Abs/ppm | | |
| 0 - 2.50 mg/l | 0.0676 | 0.0550 to 0.0800 | Pass |
| 2.50 - 5.00 mg/l | 0.0628 | 0.0600 to 0.0800 | Pass |
| 5.00 - 7.69 mg/l | 0.0625 | 0.0550 to 0.0800 | Pass |
| 7.69 - 10.0 mg/l | 0.0537 | 0.0500 to 0.0700 | Pass |

Notes

- Store between 2 & 8 °C
- *Stable on Board TON 1 - 1 day
- * TON 2 & TON 3 - 5 days
- Keep Tightly Sealed when not in use.
- Do not return unused portions to container
- Barcode ID- TON 1 103; TON 2 104; TON 3 105
- * Based on results obtained in our Research Facility

Alan Cooney
A.Cooney

An ISO 9001:2008 Registered Company

Annacotty Business Park, Annacotty, Limerick, Ireland.

Telephone 061 358190

Facsimile 061 358191

Hazard Identification

Labelling according to Regulation (EC) No 1272/2008


 Signal Word **Warning**
H315 Causes skin irritation

H319 Causes serious eye irritation

 see full MSDS on www.serosep.com

Precautionary statements

P280;P264;P321

P305+P351+P338

P302+P352; P332+P313

P362+364; P337+P313