

## Certificate of Performance

### Phosphate 2 (Phos 2) Reagent

**Product Code** KR107

**Lot Number** 4351

**Date** 13.03.2020


**Expiry Date** 12.03.2021

Reagent	Lot Number
Phos 1	4355
<b>Phos 2</b>	<b>4351</b>
<b>Test Method</b>	WC SOP 21
<b>Test Instrument</b>	Konelab 20
<b>Serial Number</b>	P4919028

Standard Value	Lot #	Abs.
0		0.001
0.5	P20 03 20	0.149
1	P20 03 20	0.302
1.5	P20 03 20	0.456
2	P20 03 20	0.595

\*serial dilution

### Calibration Curve



Test	Result	Specification	Pass / Fail
Intercept	0.002	-0.005 to 0.010	<b>Pass</b>
Sensitivity	Abs/ppm 0.299	0.260 to 0.320	<b>Pass</b>
Regression	1.000	0.990 to 1.001	<b>Pass</b>

**Notes** Store between 2 & 8 °C  
 \*Stable on Board 5 days  
 Barcode ID = 107  
 \*Based on results obtained in our Research Facility

*Alan Coonerty*  
 A. Coonerty

An ISO 9001:2008 Registered Company  
 Annacotty Business Park, Annacotty, Limerick, Ireland.  
 Telephone 061 358190 Facsimile 061 358191

### Hazard Identification

Labelling according to Regulation (EC) No 1272/2008  
 This reagent is classed as non hazardous.