

## Certificate of Performance

### Sulphate Reagent

**Product Code** KR102

**Lot Number** 4359

**Date** 13.03.2020


**Expiry Date** 24.02.2021

Reagent	Lot Number
Sulphate	4359
Test Method	WC SOP 24
Test Instrument	Konelab 20
Serial Number	P4919028

Standard Value	Lot #	Abs.
0		0.005
250	Sulp 4000 12 19	0.149
500	Sulp 4000 12 19	0.300
800	Sulp 4000 12 19	0.457
1000	Sulp 4000 12 19	0.538

\*serial Dilution

### Calibration Curve



Test	Result	Specification	Pass / Fail
Intercept	0.005	-0.005 to 0.010	Pass
Sensitivity Abs/ppm			
0 - 250 mg/l	0.0006	0.0003 to 0.0006	Pass
250 - 500 mg/l	0.0006	0.0003 to 0.0007	Pass
500 - 800 mg/l	0.0005	0.0004 to 0.0007	Pass
800 - 1000 mg/l	0.0004	0.0004 to 0.0008	Pass

**Notes** Store between 2 & 8 °C  
 \*Stable on Board 5 days  
 Barcode ID = 102  
 \* Based on results obtained in our Research Facility

*Allen Coonerty*  
 A. Coonerty

An ISO 9001:2008 Registered Company

Annacotty Business Park, Annacotty, Limerick, Ireland.

Telephone 061 358190

Facsimile 061 358191

### Hazard Identification

Labelling according to Regulation (EC) No 1272/2008

This reagent is classed as non hazardous.