

Certificate of Performance

Hard 2

Product Code KR118

Date 12.03.2020


Expiry Date 28.02.2021

Reagent	Lot Number
Hard 1	A19-159
Hard 2	4382
Hard 3	A20-01
Method	WC SOP 25
Test Instrument	Konelab 20
Serial Number	P4919028

Std. Value	Lot #	Abs.
20.00	Hard 1000 11 19	-0.115
50.00	Hard 1000 11 19	-0.300
100.00	Hard 1000 11 19	-0.626
142.00	Hard 1000 11 19	-0.846
200.00	Hard 1000 11 19	-1.062

*serial dilution

Calibration Curve



Test	Result	Specification	Pass/Fail
Sensitivity			
20 - 50 mg/l	-0.0062	-0.0050 to -0.0080	Pass
50 - 100 mg/l	-0.0065	-0.0050 to -0.0080	Pass
100 - 142 mg/l	-0.0052	-0.0040 to -0.0080	Pass
142 - 200 mg/l	-0.0037	-0.0030 to -0.0070	Pass

Notes Store between 2 & 8 °C

Keep Tightly Sealed when not in use.

Do not return unused portions to container

Barcode ID- Hard 1 117; Hard 2 118; Hard 3 119

* Based on results obtained in our Research Facility



A. Coonerty

An ISO 9001:2008 Registered Company
 Annacotty Business Park, Annacotty, Limerick, Ireland.
 Telephone 061 358190 Facsimile 061 358191

Hazard Identification

Labelling according to Regulation (EC) No 1272/2008

This reagent is classed as non hazardous.