

## Certificate of Performance

### Hard 1

Product Code KR117

Date 08.04.2020


Expiry Date 31.03.2021

Reagent	Lot Number
Hard 1	4388
Hard 2	A19-133
Hard 3	A19-134
Method	WC SOP 25
Test Instrument	Konelab 20
Serial Number	P4919028

Std. Value	Lot #	Abs.
20.00	Hard 1000 11 19	-0.105
50.00	Hard 1000 11 19	-0.297
100.00	Hard 1000 11 19	-0.640
142.00	Hard 1000 11 19	-0.882
200.00	Hard 1000 11 19	-1.152

\*serial dilution

### Calibration Curve



Test	Result	Specification	Pass/Fail
Sensitivity			
20 - 50 mg/l	-0.0064	-0.0050 to -0.0080	Pass
50 - 100 mg/l	-0.0069	-0.0050 to -0.0080	Pass
100 - 142 mg/l	-0.0058	-0.0040 to -0.0080	Pass
142 - 200 mg/l	-0.0047	-0.0030 to -0.0070	Pass

Notes Store between 2 & 8 °C

Keep Tightly Sealed when not in use.

Do not return unused portions to container

Barcode ID- Hard 1 117; Hard 2 118; Hard 3 119

\* Based on results obtained in our Research Facility



A. Coonerty

An ISO 9001:2008 Registered Company  
 Annacotty Business Park, Annacotty, Limerick, Ireland.  
 Telephone 061 358190 Facsimile 061 358191

### Hazard Identification

Labelling according to Regulation (EC) No 1272/2008

This reagent is classed as non hazardous.