

Certificate of Performance

Ammonia 1 Reagent

Product Code **KR108**
 Lot Number **4396**

Date 20.05.2020


Expiry Date 31.03.2021

Reagent	Lot Number
AMM 1	4396
Amm 2	4408
Test Method	WC SOP 22
Test Instrument	Konelab 20
Serial Number	P4919028

Std. Value	Lot #	Abs.
0.00		0.027
0.25	AMM 10 02 20	0.224
0.50	AMM 10 02 20	0.443
0.77	AMM 10 02 20	0.662
1.00	AMM 10 02 20	0.854

*serial dilution

Calibration Curve



Test	Result	Specification	Pass / Fail
Intercept	0.0236	0.000 to 0.025	Pass
Sensitivity Abs/ppm	0.8301	0.750 to 0.890	Pass
Regression	0.9998	0.990 to 1.0001	Pass

Notes
 Store between 2 & 8 °c
 *Stable on Board 1 day
 Barcode ID = 108
 *Based on results obtained in our Research Facility

Alan Coonerty
 A.Coonerty

An ISO 9001:2008 Registered Company
 Annacotty Business Park, Annacotty, Limerick, Ireland.
 Telephone 061 358190 Facsimile 061 358191

Hazard Identification

Labelling according to Regulation (EC) No 1272/2008
 This reagent is classed as non hazardous.