

Certificate of Performance

Sulphate Reagent

Product Code KR102

Lot Number 4407

Date 14.05.2020


Expiry Date 12.05.2021

Reagent	Lot Number
Sulphate	4407
Test Method	WC SOP 24
Test Instrument	Konelab 20
Serial Number	P4919028

Standard Value	Lot #	Abs.
0		0.003
250	Sulp 4000 12 19	0.142
500	Sulp 4000 12 19	0.298
800	Sulp 4000 12 19	0.463
1000	Sulp 4000 12 19	0.554

*serial Dilution

Calibration Curve



Test	Result	Specification	Pass / Fail
Intercept	0.003	-0.005 to 0.010	Pass
Sensitivity Abs/ppm			
0 - 250 mg/l	0.0006	0.0003 to 0.0006	Pass
250 - 500 mg/l	0.0006	0.0003 to 0.0007	Pass
500 - 800 mg/l	0.0006	0.0004 to 0.0007	Pass
800 - 1000 mg/l	0.0005	0.0004 to 0.0008	Pass

Notes Store between 2 & 8 °C
 *Stable on Board 5 days
 Barcode ID = 102
 * Based on results obtained in our Research Facility

Allen Coonerty
 A. Coonerty

An ISO 9001:2008 Registered Company

Annacotty Business Park, Annacotty, Limerick, Ireland.

Telephone 061 358190

Facsimile 061 358191

Hazard Identification

Labelling according to Regulation (EC) No 1272/2008

This reagent is classed as non hazardous.