

## Certificate of Performance

### Hard 3

Product Code KR119

**Date** 02.06.2020


**Expiry Date** 02.06.2021

Reagent	Lot Number
Hard 1	A20-68
Hard 2	A20-79
Hard 3	<b>4416</b>
Method	WC SOP 25
Test Instrument	Konelab 20
Serial Number	P4919028

Std. Value	Lot #	Abs.
20.00	Hard 1000 11 19	-0.116
50.00	Hard 1000 11 19	-0.313
100.00	Hard 1000 11 19	-0.661
142.00	Hard 1000 11 19	-0.947
200.00	Hard 1000 11 19	-1.204

\*serial dilution

### Calibration Curve



Test	Result	Specification	Pass/Fail
Sensitivity			
20 - 50 mg/l	-0.0066	-0.0050 to -0.0080	Pass
50 - 100 mg/l	-0.0070	-0.0050 to -0.0080	Pass
100 - 142 mg/l	-0.0068	-0.0040 to -0.0080	Pass
142 - 200 mg/l	-0.0044	-0.0030 to -0.0070	Pass


**Notes** Store between 2 & 8 °C

Keep Tightly Sealed when not in use.

Do not return unused portions to container

Barcode ID- Hard 1 117; Hard 2 118; Hard 3 119

\* Based on results obtained in our Research Facility

  
 A. Coonerty

An ISO 9001:2008 Registered Company  
 Annacotty Business Park, Annacotty, Limerick, Ireland.  
 Telephone 061 358190 Facsimile 061 358191

### Hazard Identification

Labelling according to Regulation (EC) No 1272/2008

This reagent is classed as non hazardous.